Get to know MedMinder







Falls and incorrect use of medication often force seniors out of their homes and into institutionalized care.

MedMinder offers wireless devices that provide an all-in-one medication management and an emergency alert system.

This around-the-clock medication reminder and monitoring system is simple, convenient, and reliable.

Our emergency alert system is embedded directly into the device and communicates with pendants worn as a necklace or watch around the home. With the push of a button, our trained and certified monitoring professionals are ready to help when it is needed most.




Designed with ease of use in mind, our simple yet sophisticated devices are painless to program and user friendly.

MedMinder pill dispensers mirror pharmacy pill organizers, this creates an immediate familiarity. The cellular modem embedded in the box allows seamless programming and monitoring to support the patient and notify all interested parties.

Patients who may be resistant to new technology will feel comfortable using a MedMinder pill dispenser. There are no digital read-outs, buttons or complicated programming required for the MedMinder. There is no need for a computer, any form of Internet access or technical skills.





Reminders

Users are prompted to take or refill their medications through a series of customizable visual and auditory alerts. Set up the reminders to blink, beep or play a custom audio message. If the user doesn't respond to these reminders they will receive a phone call. If there is no response the family or caregiver are notified by phone, email or text message ensuring that everyone stays informed.

Remote Progamming And Monitoring

MedMinder is equipped with a self-contained cellular modem that links to our central server. Using our secured website users and caregivers can easily program their unit's medication schedule remotely, set up preferred notifications and review reports.

Medical Alert

With the push of a button the pill dispenser will open a two-way voice channel with the professionals at the UL certified monitoring center.

The wireless pendants are waterproof, feature 800 feet of transmission range, a long-life battery, and have a necklace and a male or female watch band design.

Bluetooth

MedMinder devices are equipped with bluetooth connectivity and can serve as gateways for other telehleath devices.



REMOTE SUPPORT SERVICES

www.RSSMonitoring.com

1-888-325-2094

(Re)filling

MedMinder offers two easy options to fill the medication compartments:

- 1. **Manually:** Simply lifting the lockable main lid gives access to all 28 compartments.
- 2. **Pre-fill tray:** MedMinder supplied pre-fill trays can be filled by a pharmacist, caregiver or family member. Once filled they can be placed directly into MedMinder.

Locking

The Jon pill dispenser is for those at greater risk for taking the wrong medication at the wrong time. All compartments are locked and only the specified compartment will unlock at the right time. The Jon features an "un-locking / re-locking schedule" where the dosage compartments can be setup to un-lock or re-lock before or after the scheduled dosage to allow more flexibility for the user.

Custom Greetings

Caregivers can record voice messages at the MedMinder portal and our devices will play them so the user will hear a familiar voice. Custom greetings can be set as medication reminders, positive feedback, or even one-time personalized messages for special occasions. Caregivers can also type a personalized message that MedMinder converts to voice format.



Reporting

MedMinder data is designed for future integration with existing medical records. Integrating MedMinder's recorded medication activity with existing records will allow better data correlation and improved quality of care.

Feature	Maya	Jon
Simple Design	Ø	<
Flexible Programming	Ø	\checkmark
Numerous Reminders	\bigcirc	$ \bigcirc $
Multiple Caregivers	Ø	Ø
Large Capacity	Ø	igstar
Multi Weeks	Ø	\checkmark
Cellular Connection	⋖	$ \checkmark $
Pre-Fill Tray Option	\bigcirc	\checkmark
Locking Mechanism		Ø
Medical Alert		