

Patient Instructions: How To Attach an Applicator to a Cream Tube & Finding the Right Dose

Step 1:



Remove the protective seal from the tube of cream. (This step is only required before first use)

Step 2:



Attach the applicator to the tube by **firmly pressing** the tip onto the opening while **twisting with steady force** until it is securely in place.

Step 3:



Once firmly attached, slowly squeeze the tube with an upward motion in order to fill the applicator. Continue until the plunger reaches the prescribed dose marking. The line should align with the top edge of the applicator.



Step 4:



Unscrew the filled applicator from the tube. It is now ready to use according to your prescription directions.

Step 5:



After use, replace the cap on the tube. Wash the applicator with warm water, allow it to air dry, and store it for your next dose.