

An in-home test designed to get accurate results.

Not only is the iHealth COVID-19 Antigen Rapid Test easy and non-invasive, it's a small engineering beauty. The test has a highly-absorbent swab that is individually sealed to help prevent contamination so you can view the results of your in-home test with confidence.




Highly-Absorbent swab

Made of PU foam, the absorbent tip of the swab makes sample collection more efficient than a regular swab.



Double-Cap Security Design

The double-cap design allows you to safely store the liquid and then precisely drip the mixed sample onto the test card.




Easy-to-Understand Results

The unique test card shortens the time to get your results and makes them easy to read. If only the C line shows up, the test is negative. If both C and T lines show up, the test is positive.

Independently-sealed packaging ensures safe and hygienic testing.

The swabs and test cards of iHealth COVID-19 Antigen Rapid Test are sealed and independently packaged to ensure the safety and hygiene of the product.



How to use the iHealth COVID-19 Antigen Rapid Test:

