

HRTICKER®

"THE CURE FOR THE COMMON APPLICATOR"

INSTRUCTIONS FOR PATIENTS

Please read all information carefully before using the HRTicker.

STEP 1



Remove the safety lid by placing your thumb under the thumb tab and pushing upward.

PRIMING THE UNIT

If you begin to turn the applicator and the cream does not come out, it may be an indication that you need to prime the device. **Priming:** turn the applicator base clockwise until the cream begins to dispense out. (It may take several turns.) After the cream comes out, use a tissue to wipe off the cream. Your device is now primed and ready to use!

STEP 2



Turn the applicator base clockwise to apply your specific dose, as instructed by your pharmacist or doctor. Use as many clicks as instructed.

STEP 3



Always follow your doctor or pharmacist's instructions on how to apply the medication. If desired, apply the cream directly with the HRTicker by rubbing the dispenser cap against the area of skin as shown.

STEP 4



Prevent evaporation of the next dose by completely securing the safety lid. (TIPS: Some of the water that makes the bulk of the cream may evaporate without the safety lid.)

STEP 5



When the top part of the plunger reaches the extruded word "REFILL" (as shown above), it is time to refill your medication. There are approximately 6-7 grams remaining.

Universal Arts Compounding Pharmacy

14350 NW 56 Court, #114 Miami, FL 33054 Tel: 305-556-2673

Fax: 305-556-9749 Website: www.uaprx.com Email: compoundrx@uaprx.com