

AF3

LARGE WIPE

Sani-Cloth®

GERMICIDAL DISPOSABLE WIPE



PRODUCT DESCRIPTION

Effective against 45 microorganism in 3 minutes.

Protecting your staff and patients has never been easier! Alcohol-and fragrance-free formula is ideal for use around patients, residents and staff, especially those with respiratory sensitivities.

Disinfection Level

- EPA registered intermediate level disinfectant

Applicator

- Thick, strong wipe has high fluid capacity and is pliable to get into crevices and corners of equipment
- Pre-moistened wipes are ready to use

Overall Contact Time

- 3 minutes

Efficacy

- Bactericidal, Tuberculocidal, Virucidal

Efficacy Claims

- 45 microorganisms, including multi-drug resistant organisms

Formulation

- Alcohol free quaternary ammonium

Formulation does not contain

- Fragrance, alcohol, bleach, acid, phenol, acetone and ammonia

Packing features

- Deep well lid provides ample space to store and access next wipe
- Prominent contact time for easy identification and improved compliance with accreditation standards
- Categorization of pathogens allows for quick identification of relevant efficacy claims
- Icon clearly communicates key safety information on both label and lid

Compatibility

- Compatible with a broad range of surfaces and equipment in healthcare

Healthcare Use Sites

- May be used on hard nonporous surfaces of: bed railings, blood glucose meters (glucometers) cabinets, carts, chairs counters, dental unit instruments trays, exam tables, gurneys, isolettes IV poles, stethoscopes, stretchers, tables, telephones, toilet seats, diagnostic equipment, patient monitoring equipment, patient support and delivery equipment and MANY MORE! Please refer to the product / master label to get the full listing of use sites.

Compliance

- Meets CDC, OSHA and Infection Control Guidelines

One Step

- One step cleaner and disinfectant

Tools and Accessories

- Point of care accessories and tools are available to help improve staff compliance to facility protocols

EPA LIST N #9480-9

Packing.....	160 Count Wipes
Material	Non-Woven Fabric Rolls
Fabric	Embossed Spunlace
Weighs	40 gsm
Wipe Size	6.0" x 6.75"
KILL TIME	3 Minutes

ACTIVE INGREDIENTS

Alkyl: (68% C12, 32% C14)	
Dimethyl Benzyl Ammonium Cl:	0.14%
Alkyl: (60% C14.1, 30% C14.2, 5% C12, 5% C18)	
Dimethyl Benzyl Ammonium Cl:	0.14%
OTHER INGREDIENTS	99.72%
TOTAL.....	100.00%

PACKING

Units / Cases	12
Cases / Pallets	72 (864)
Pallets / 40 HQ Container	20 (24000)
Unit Weight.....	0.75 kg
Case Gross Weight	1.0 lb ; 11.04oz
MOQ.....	24,000 Units

CONTACT

sales@unityhealth.us



UNITY
HEALTH