



**True Compounding Pharmacy**  
634 C Pine Ridge Dr, West Columbia, SC 29172  
Phone: 803-888-2990  
Hours: Mon–Fri: 9 a.m.–5 p.m. | Sat–Sun: Closed

## The Topi-Click Device



The Topi-Click™ is a user-friendly, mess-free device that helps you apply your prescribed cream accurately and consistently. Each "click" dispenses a measured dose of your medication — no guessing needed!

### Benefits:

- ✓ Easy to use — just click and apply
- ✓ Accurate dosing every time
- ✓ Less mess, no wasted medication
- ✓ Travel-friendly and discreet

## How to Use Your Topi-Click™ Applicator:

1. Before and after applying your medication, wash your hands with soap and water.
2. Lift the cap off the Topi-Click Device by placing your thumb under the tab on the outer cap.
3. After opening the Topi-Click, wipe off the applicator dome to ensure you are getting a consistent dose.
4. Ensure that the indicator lines on the base and the tube are aligned. Then twist the base to the number of clicks specified on your prescription label. After each click, the indicator lines should align.
5. Use the dome top of the device to rub the cream directly onto clean, dry skin.
  - Common application areas include inner thighs, inner arms and abdomen. Application sites should be rotated, unless otherwise instructed.
  - The application area should be about the size of your palm or a smaller dinner plate. Rub using the dome top of the device for 2 minutes until completely absorbed.
6. Wipe off any residual medicine remaining on the dome to ensure you are getting a consistent dose and prevent any product from getting stuck on the cap. Replace cap.



## Important Reminders:

- **Cream can be transferred to other people, children or pets – choose a site that is covered and wash hands after applying.**
- Keep out of reach of children and pets.
- Store at room temperature, away from heat or sunlight.
- Do not share your cream with others.
- Do not apply other creams or moisturizers to the same area – as they can interfere with absorption.
- Once the bottom of the plunger reaches the line below “Refill” there are about 32 clicks remaining and once it reaches the line on top of the word “Refill” there are about 12 clicks remaining.