



QUIDEL®

Sofia[®]2

Accurate, objective, automated results you can rely on.

- Accurate:** Excellent sensitivity and specificity compared to molecular and viral culture methods.
- Objective:** Fluorescent technology guarantees an objective result with no potential for subjective interpretation.
- Easy:** Simple workflows with as few as one button press and ~1 minute of hands-on time.
- Fast:** Results in 3-15 minutes (differs by assay).
- Convenient:** Small instrument size and room-temp storage of kits are suitable for any office or lab.
- Efficient:** Automatically stores patient results – results never expire! View results from your computer using Virena.



WALK AWAY mode lets the Test Cassette incubate inside Sofia 2 and reports a result automatically upon completion.

READ NOW mode lets you incubate multiple Test Cassettes on the bench top and quickly analyze them one by one in Sofia 2.

Choice without compromise.

