



There's a CONTOUR®NEXT Meter for Everyone

Recommend the CONTOUR®NEXT GEN

A great choice for patients who prefer a highly accurate,^{1,2,3} easy-to-read and comfortable-to-handle meter!

Easy to Use

- Large display screen and rubber buttons
- Readings in 5 seconds
- Second-Chance® sampling allows your patients up to 60 seconds to reapply more blood if needed
- Target range indicator light immediately identifies blood glucose readings as above, below or within target range



Remarkable Accuracy^{1,2,3}

- The CONTOUR®NEXT GEN system has been shown to deliver highly accurate blood glucose readings



Seamlessly Connects to the Free CONTOUR®DIABETES App

- Reports can be shared with patient's healthcare professional



It's All About Choice!

All CONTOUR®NEXT meters use CONTOUR®NEXT test strips



References: 1. CONTOUR®NEXT GEN BGMS User Guide, Rev 9/20. 2. Christiansen, M. P. (2017). A New, Wireless-enabled Blood Glucose Monitoring System That Links to a Smart Mobile Device: Accuracy and User Performance Evaluation. *Journal of Diabetes Science and Technology*, 11(3), 567-573. DOI: 10.1177/1932296817691301. 3. Bernstein, R. et al. (2013). A New Test Strip Technology Platform for Self-Monitoring of Blood Glucose. *Journal of Diabetes Science and Technology*, 7(5), 1386-1399.