


POLOXAMER

Poloxamer is a great way to administer medication to animal ear canals. It utilizes Reverse Thermodynamics, which means it is a liquid when cold and a semi-solid when warm.



POLOXAMER VET EAR PLUGS

- Used for otitis externa, typically in dogs
- Active pharmaceutical ingredients, such as antibiotics, antifungals, steroids, or combinations of these, are mixed into the poloxamer base
- Pre-packed in 5 mL syringes with applicator
- Applied once in-clinic and left in place for 5-7 days
- Water soluble – washes out with water irrigation



ACTIVE PHARMACEUTICAL INGREDIENT EXAMPLES

Enrofloxacin 2%
Ketoconazole 2%
Gentamicin 0.3%
Triamcinolone 0.25-0.5%

Betamethasone Valerate 0.12%
Clotrimazole 1%
Combination of the above

For more information contact your nearest NuCara Compounding Pharmacy:

Ames, IA:

Ph: 515.292.3604

Coralville, IA:

Ph: 319.354.6006

Ottumwa, IA:

Ph: 641.683.8168

Pleasant Hill, IA:

Ph: 515.264.1503



Waterloo, IA:

Ph: 319.236.8891

Elmhurst, IL:

Ph: 630.834.9768

Austin, TX:

Ph: 512.454.9923



www.nucara.com